

letto bilancia mod.
bed-scale mod.

OMEGA T

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The bed-scale is supplied with: built-in weighing system; one pushbutton panel for the patient and one mobile console for the caregiver for moving the backrest, double translating pelvic section, variable height with columns, trendelenburg position, counter-trendelenburg position, cardiological seat position and more upright therapy chair position, with resetting of all movements; rapid backrest and RCP release with pedals; ABS bed base and parts; built-in extension; antimicrobial treatment; carbon steel structure subjected to mechanical and chemical treatments; wheels complete with braking system.

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Bed fitted with weighing system, with bed base in 3 jointed sections and 1 fixed section, suitable for cardiac massage; made with easily removable ABS panels. The particular shape of the structure supporting the bed guarantees a large radio-transparent zone. The bed has a pushbutton panel for the patient and a command console exclusively for staff use. All the commands on the patient's pushbutton panel may be enabled or disabled with the console. Buttons which are enabled are identified by led lights. The pushbutton panel also has a button for simultaneously moving the backrest and auto-contour thigh-rest.

The backrest and thigh rest automatically translate during rotation. Correct posture not only facilitates breathing but reduces the possibility of bed sore formation. Sophisticated electronic control of bed movements makes it possible to work the patient gradually into the 'lounge chair' position. There is also a more vertical chair position, known as the 'therapy chair', which facilitates ventilation of the patient's lungs.

The bed has two pedals, positioned on either side, for rapidly releasing the backrest in an emergency; when the pedal is pressed down it mechanically releases the backrest (RCP), levels the bed base surface and positions it at minimum height to permit cardiac massage, while if lifted up it electrically activates the Trendelenburg position (-18°) while automatically levelling the surfaces. The pedals on the base near the head area allow staff to work on the patient with free hands.

The Trendelenburg and counter-trendelenburg movements may be electrically controlled with the staff-only control console. The leg rest section has a built-in extension.

The headboard is designed to permit use as a handle for pushing the bed; the headboard and footboard are removable; at all four corners of the bed there are housings for insertion of a beam for lifting the patient or an IV support rod and plastic bumpers.

The braking system has 4 pedals (one at each corner of the bed), each of which locks and releases the bed's wheels at the same time. An acoustic alarm will sound if the bed is plugged in without the brakes on the wheels. To facilitate movement along longer routes, the direction of the front wheel may be locked. The lower central part of the base is completely free to permit insertion of a beam for lifting the patient.

The innovative four section overturning sides, with shock absorbed vertical travel, have a translation synchronisation system that follows the motion of the bed base surface to permit optimal containment of the patient in all bed positions.

The bed structure is made with carbon steel subjected to special mechanical and chemical treatments to guarantee the best possible aesthetic and functional results and durability, guaranteeing ongoing high resistance to oxide formation and complete absence of rust. The structure is subjected to a certified antimicrobial treatment to eliminate the proliferation of bacteria. The total absence of corners and gaps makes it easy and safe to use and keeps it clean at all times.

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